



sartorius stedim  
biotech

## Automated Glucose Control



Application  
Note

#07

#08

#09

#10

#11

turning science into solutions

Authors:  
Stuart Tindal, Sebastian Ruhl, Diana Hesse

### Introduction



One of the goals of mammalian upstream bioprocess development is to handover a process to clinical manufacturing that is robust, safe and in the end produces enough material to meet the studies demand. These process development teams are increasingly being asked to do more and in less time. Thus, automated feeding strategies are becoming more and more popular as a way of reducing the development efforts and as a side benefit have shown to improve process performance and product attribute consistency.

Cell cultivation's primary feed component is glucose which is used as a starting point for all growth and energy pathways within the cells. Ensuring that the cells do not have too much or too little, enables them to grow fast and maximize the product secretion. Additionally, it has been noted that an excess in glucose concentration has an effect on the glycosylation rate of secreted soluble proteins. Thus, controlling the available glucose would improve the consistency of glycosylation patterns on final products and may improve the quality.

Within this application note a stepwise example method on how to establish glucose feed control is documented. The materials and methods show how this was done using Sartorius Stedim Biotech hardware and equipment. An overview of the results and discussion of the data is given in order to highlight some of the key benefits of applying this method to other processes.

### Material & Methods

A BIOSTAT® B-DCU was used as bioreactor control system. The BIOSTAT® B-DCU is a bioreactor for advanced process optimization and characterization featuring the option of a glucose concentration controller. The integration of the BioPAT® Trace allows a real-time monitoring of the glucose level enabling a user defined glucose setpoint within the BIOSTAT® B-DCU control software. The first two runs, insitu glucose concentration was maintained by discontinuous bolus feeds. Subsequently the third run utilized a glucose setpoint controller using defined PID settings and an internal, speed-controlled peristaltic pump.





The BioPAT® Trace was set up for dialysis mode glucose & lactate measurement at a 20 minutes sampling rate. It was set to a fully automated and self-calibrating protocol as part of the schedule tab function. Each 17 day run used 12 L of transport buffer, two calibration solutions (high: 10 g/L glucose; low: 1 g/L glucose) and a 10 L waste container attached to the tube set. The dialysis probe was prepared, filled and installed into the UniVessel® Glass 5 L prior to sterilization and connected to the BioPAT® Trace tube set and primed for analysis before inoculation.

The UniVessel® Glass 5 L was equipped with two 3-blade segment impellers for low shear stress and good homogenization of the cell broth. The blade angle was 30° and set to down pumping. A ring sparger with holes faced up was used for all trials. Additionally, the vessel was equipped with several ports for feeding, a classical pH sensor, pO<sub>2</sub> sensor, dialysis probe, exhaust cooler and gas filters.

To evaluate the BioPAT® Trace integration in the BIostat® B-DCU a CHO fed-batch process was used. The 17 day cultivation comprises of a 3 day batch phase and a 14 day fed-batch phase. After the inoculation with  $0.3 \times 10^6$  cells/mL the peak viable cell density (VCD) is typically reached at day 8 with  $25 - 30 \times 10^6$  cells/mL and a viability of 99 %. After the following 9 day dying phase the VCD should be above  $10 \times 10^6$  cells/mL with a viability of more than 50 % at the point of harvest.

The bolus feeding from day 3 comprises feed medium A (FMA), feed medium B (FMB) and a highly concentrated glucose solution (400 g/L). The fed amount of FMA and FMB is constant throughout the complete fed-batch phase. Typically on day 7, additional glucose is needed to maintain a glucose concentration of at least 3 g/L in the cell broth. The feeding process is automated using balances and pumps connected to the digital control unit (DCU) and S88 recipe in the SCADA software BioPAT® MFCS.

Analytics were an essential part of this evaluation. Among other things, offline glucose and lactate measurements were performed with the Radiometer ABL 800 basic. VCD and viability were analyzed with a Cedex HiRes.

A total of three different trials were performed using the following characteristics:

1. Monitoring  
Trace technology was solely used to monitor the online glucose trend throughout the cultivation.
2. Online value replaces offline measurement  
During the CHO culture a daily sample is taken. The offline glucose value is used to fill to a certain glucose concentration in the bioreactor. In trial 2 the glucose online value is used instead of the offline value for this feeding procedure.
3. Glucose control  
During the third trial glucose control at 6 g/L was activated after 6 days. The usual bolus feed of FMA was modified to continuous feeding to reduce daily peaks in the glucose concentration. After glucose reached a concentration below 6 g/L, a controlled feed of glucose solution (400 g/L) was initiated to maintain the glucose concentration at 6 g/L.

Feeding scheme for trial 1 and 2:




Figure 1: Feeding scheme 1 and 2

For trials 1 and 2 the feeding scheme was identical. During the batch phase no medium was fed to the culture. On day 3 a daily bolus feed of FMA and FMB started. After glucose fell below 6 g/L an additional feed of glucose (400 g/L) started.




Figure 2: Feeding scheme 3

In trial 3 the bolus feeding of FMA was modified to continuous feeding to remove daily peaks in the glucose concentration. After glucose reached a concentration below 6 g/L, a controlled feed of glucose solution (400 g/L) was initiated to maintain the glucose concentration at 6 g/L.

## Results & Discussion

The batches ran sequentially on the same BIOSTAT® B-DCU and BioPAT® Trace systems. The results are shown in figure 3 and 4 to demonstrate the reproducibility and conformance of the process batches (compared to the historical golden batch). The glucose monitoring of the BioPAT® Trace and subsequent glucose controller are shown in figures 5-7.




Figure 3: VCD trend trial 1-3

Trial 1 and 2 fit the golden batch trend well. Due to a general process modification (continuous FMA feed) the VCD trend of trial 3 has a slightly reduced peak VCD. Final VCD and viability fits well or exceeds golden batch values for all trials.




Figure 4: IgG trend

Comparable productivity in all three trials was demonstrated by daily measurements of the product concentration (IgG concentration, figure 4). Due to a lower cell growth in trial three a slightly lower product concentration was achieved in this run.




Figure 5: Glucose and lactate trend trial 1

Despite of a small offset in the first 24 hours the online and offline measurements of glucose and lactate fit well during the cultivation.




Figure 6: Glucose and lactate trend trial 2

For trial 2 the online and offline glucose and lactate measurements showed good comparability.






Figure 7: Glucose and lactate trend trial 3

Also in trial 3 the online and offline trends fit well for glucose and lactate. At day 13 and 15 an unnecessary recalibration was performed (operational mistake).

The transition from glucose monitoring to glucose control requires a change of the process. The effects of tighter glucose control are known to be beneficial to the process performance and product quality (reference An Zhang et al). Considering this, changing a process using non tried and tested technology is undesired. Thus, learning the system and building confidence around the measurement output (comparing to offline sampling) and ensuring the device's robustness is generally needed before making that change. The first two process runs show the high resolution monitoring capabilities of the integrated BIOSTAT® B-DCU and BioPAT® Trace whereas the third run engages the glucose controller. The overall user operational interactions with the bioreactor related to sampling glucose are removed. This means the normal hours and out of normal working hours needed to maintain the BIOSTAT® B-DCU system in a glucose controlled state are reduced and more manageable. In table 1 a summary of the time savings are presented, in addition all user interaction work with the system was done within the initial stages of the BIOSTAT® B-DCU setup.

Table 1: FTE resource balance per batch when using the BioPAT® Trace for glucose

Process Step	Hours
Setup of trace probe and electronics	2
Saved time for sampling	9 (18 × 0.5)
Saved time	7

### Summary & Conclusion

The stepwise integration of the glucose control technology to a CHO fed-batch process was successful and maintained an inprocess steady-state glucose concentration of 6 (+/- 0.25) g/L.

The direct integration of the BioPAT® Trace technology into the BIOSTAT® B-DCU is fully functional and the utilization was user friendly, easy to establish and dynamically variable when required. As the system package is coming from Sartorius Stedim Biotech completely, the lifecycle of the glucose controller is managed by us through all scales of development and commercial manufacturing. This ensures the benefits uncovered in development can be transferred to commercial production for the complete product lifecycle.

### Outlook

Further trials with varying, automatically controlled glucose concentrations will be performed in order to improve product concentration and product quality attributes. Additionally, the operator time saving (7 hours per batch) and manageable working hours makes steps towards a walk away bioreactor.

## Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

### Europe

#### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5-7  
34302 Guxhagen

Phone +49.5665.407.0  
Fax +49.5665.407.2200

#### France

Sartorius Stedim FMT S.A.S.  
ZI des Paluds  
Avenue de Jouques - CS 91051  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.845619

Sartorius Stedim France SAS  
ZI des Paluds  
Avenue de Jouques - CS 71058  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.846545

#### Austria

Sartorius Stedim Austria GmbH  
Modecenterstrasse 22  
1030 Vienna

Phone +43.1.7965763.18  
Fax +43.1.796576344

#### Belgium

Sartorius Stedim Belgium N.V.  
Rue Colonel Bourg 105  
1030 Bruxelles

Phone +32.2.756.06.80  
Fax +32.2.756.06.81

#### Hungary

Sartorius Stedim Hungária Kft.  
Kagyló u. 5  
2092 Budakeszi

Phone +36.23.457.227  
Fax +36.23.457.147

#### Italy

Sartorius Stedim Italy S.r.l.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41  
Fax +39.055.63.40.526

#### Netherlands

Sartorius Stedim Netherlands B.V.

Phone +31.30.60.25.080  
Fax +31.30.60.25.099

[filtratie.nederland@sartorius-stedim.com](mailto:filtratie.nederland@sartorius-stedim.com)

#### Poland

Sartorius Stedim Poland Sp. z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn

Phone +48.61.647.38.40  
Fax +48.61.879.25.04

#### Russian Federation

LLC "Sartorius Stedim RUS"  
Uralskaya str. 4, Lit. B  
199155 St. Petersburg

Phone +7.812.327.53.27  
Fax +7.812.327.53.23

#### Spain

Sartorius Stedim Spain, S.A.U.  
Avda. de la Industria, 32  
Edificio PAYMA  
28108 Alcobendas (Madrid)

Phone +34.913.586.098  
Fax +34.913.589.623

#### Switzerland

Sartorius Stedim Switzerland AG  
Ringstrasse 24 a  
8317 Tagelswangen

Phone +41.52.354.36.36  
Fax +41.52.354.36.46

#### U.K.

Sartorius Stedim UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737159  
Fax +44.1372.726171

#### Ukraine

LLC "Sartorius Stedim RUS"  
Post Box 440 "B"  
01001 Kiev, Ukraine

Phone +380.44.411.4918  
Fax +380.50.623.3162

### Americas

#### USA

Sartorius Stedim North America Inc.  
5 Orville Drive, Suite 200  
Bohemia, NY 11716

Toll-Free +1.800.368.7178  
Fax +1.631.254.4253

#### Argentina

Sartorius Argentina S.A.  
Int. A. Ávalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505  
Fax +54.11.4762.2333

#### Brazil

Sartorius do Brasil Ltda  
Avenida Senador Vergueiro 2962  
São Bernardo do Campo  
CEP 09600-000 - SP - Brasil

Phone +55.11.4362.8900  
Fax +55.11.4362.8901

#### Mexico

Sartorius de México, S.A. de C.V.  
Libramiento Norte de Tepetzotlan s/n,  
Colonia Barrio Tlacateco,  
Municipio de Tepetzotlan,  
Estado de México,  
C.P. 54605

Phone +52.55.5562.1102  
Fax +52.55.5562.2942

[leadsmex@sartorius.com](mailto:leadsmex@sartorius.com)

#### Peru

Sartorius Peru S.A.C.  
Av. Emilio Cavenecia 264 San Isidro  
15073 Lima, Perú

Phone +51.1.441 0158  
Fax +51.1.422 6100

### Asia | Pacific

#### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800  
Fax +61.3.8762.1828

#### China

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
No. 33 Yu'an Road  
Airport Industrial Park Zone B  
Shunyi District, Beijing 101300

Phone +86.10.80426516  
Fax +86.10.80426580

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
3rd Floor, North Wing, Tower 1  
No. 4560 Jinke Road  
Zhangjiang Hi-Tech Park  
Pudong District  
Shanghai 201210, P.R. China

Phone +86.21.6878.2300  
Fax +86.21.6878.2882

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
Guangzhou Representative Office  
Unit K, Building 23  
Huihua Commerce & Trade Building  
No. 80 Xianlie Middle Road  
Guangzhou 510070

Phone +86.20.37618687 | 37618651  
Fax +86.20.37619051

#### India

Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra  
Nelamangala Tq  
562 123 Bangalore, India

Phone +91.80.4350.5250  
Fax +91.80.4350.5253

#### Japan

Sartorius Stedim Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300  
Fax +81.3.4331.4301

#### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

#### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
1 Science Park Road,  
The Capricorn, #05-08A,  
Singapore Science Park II  
Singapore 117528

Phone +65.6872.3966  
Fax +65.6778.2494

#### South Korea

Sartorius Korea Biotech Co., Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, Bundang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700  
Fax +82.31.622.5799



◀ [www.sartorius-stedim.com](http://www.sartorius-stedim.com)